



# UroWIRE<sup>®</sup>

3-in-one Guidewire



# UroWIRE

## 3-in-one Guidewire

---


### With Kink-Resistant Nitinol Core

Simplified Access

Superior Handling

### Ultimate Performance

The UroWIRE guidewire facilitates access, allows for easy manipulation and protects delicate ureteroscopes. The unique design eliminates costly and time-consuming guidewire exchanges. Ideally suited for ureteroscopy, the UroWIRE guidewire eases procedures, reducing O.R. time and procedural costs.



#### Superior Handling


Super stiff midsection to facilitate instrument passage, PTFE coating reduces risk of glide out

#### Easier Access

Hydrophilic coated tip navigates tortuous anatomy and obstructed ureters

#### Scope Protection

Flexible distal end facilitates passage of delicate ureteroscopes, protecting them from damage



## UroWIRE 3-in-one Guidewire

Reorder No.	Size	Length	Qty/Box
B8634	.035"	150cm	6/BOX

Visit [www.appliedmedical.com/endourologic](http://www.appliedmedical.com/endourologic)  
for more information

© 2018 Applied Medical Resources Corporation. All rights reserved.  
Applied Medical, the Applied Medical logo design and marks designated with a © are trademarks  
of Applied Medical Resources Corporation, registered in one or more of the following countries:  
Australia, Canada, Japan, South Korea, the United States and/or the European Union.  
SC01073J

