

**With Kink-Resistant  
Nitinol Core**

- ▶ Simplified Access
- ▶ Superior Handling

# UroWIRE®

3-in-one Guidewire

**Applied**  
Medical



*A New Generation Medical Device Company*

# UroWIRE® 3-in-one Guidewire

## Easier Access

*Hydrophilic coated tip navigates tortuous anatomy and obstructed ureters.*

## Scope Protection

*Flexible distal end facilitates passage of delicate ureteroscopes, protecting them from damage.*

## Superior Handling

*Super stiff midsection to facilitate instrument passage, PTFE coating reduces risk of glide out.*

## Ultimate Performance

The UroWIRE guidewire facilitates access, allows for easy manipulation, and protects delicate ureteroscopes. The unique design eliminates costly and time-consuming guidewire exchanges. Ideally suited for ureteroscopy, the UroWIRE guidewire eases procedures, reducing O.R. time and procedural costs.

REORDER NO.	SIZE	LENGTH	QUANTITY
B8634	.035"	150cm	6/BOX



© 2013 Applied Medical Resources Corporation. All rights reserved.  
Applied Medical, the Applied Medical logo design, ProTerra, the ProTerra logo design and marks designated with a ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, the United States and/or the European Union. SC01073H

<http://contact.appliedmedical.com/>

For more information [www.appliedmedical.com/accessories](http://www.appliedmedical.com/accessories)

