




**Epix<sup>®</sup>**

Laparoscopic Instrumentation

# RELIABILITY, PRECISION AND CONTROL



## Epix Scissors

High performance blades sustain sharpness throughout duration of procedure.

Low profile design eliminates protrusions, enabling full insulation down to the blades.



## Epix Graspers

Clinically proven *latis*<sup>®</sup> pads gently engage tissue, reducing the pressure necessary to maintain traction.

Handle design incorporates the ability to grasp and lock tissue as desired, ensuring precise control.



## Epix Dissector

Curved jaws enhance the grasping, mobilizing and cauterizing of tissue.

Tapered, serrated jaws provide traction to engage tissue.



## Disposable Scissors

Reorder No.	Size	Qty/Box
-------------	------	---------

*Laparoscopic*

CB030	5mm x 35cm	10/Box
CB040	5mm x 45cm	10/Box



## Disposable Graspers

Reorder No.	Size	Qty/Box
-------------	------	---------

*Laparoscopic*

C4130	5mm x 35cm	10/Box
C4140	5mm x 45cm	10/Box



## Disposable Dissectors

Reorder No.	Size	Qty/Box
-------------	------	---------

*Laparoscopic*

CY010	5mm x 35cm	10/Box
CY020	5mm x 45cm	10/Box

Visit [www.appliedmedical.com/epix](http://www.appliedmedical.com/epix)  
for more information