



Epix[®]

Laparoscopic Suction/
Irrigation System

Epix

Laparoscopic Suction/Irrigation System

THE ULTIMATE IN COMFORT AND CONTROL

Reusable Probe Attachments



C7205



C7202/C7211



C7204

Probe Attachments

Connect and disconnect quickly and easily

Available in a variety of reusable probe lengths, diameters and tip offerings that enable a wide range of flow rates to be achieved with a single handpiece

Soft, Flexible Tubing

Drapes easily, minimizes coiling and improves control


Built-In Smoke Evacuation

Maintains visibility and eliminates the need for a separate smoke evacuation device

Color-Coded Buttons with Soft Springs

Depress smoothly and maximize user comfort

Differentiate the suction and irrigation buttons



Clear, Ergonomic Handpiece

Provides a comfortable fit and clear unobstructed view of the working channels



Epix Disposable Tubing Sets

Reorder No.	Description	Diameter	Working Length	Qty/Box
C6001	Probe & valve with dual bag tubing	5mm	33cm	12/Box
C6003	Valve only with dual bag tubing*	N/A	N/A	12/Box
C6005	Valve only with 12" tubing**	N/A	N/A	12/Box



Epix Reusable Probes

Reorder No.	Description	Diameter	Working Length	Qty/Box
C7202	Probes with irrigation holes	5mm	33cm	1/Box
C7211	Probes with irrigation holes	5mm	45cm	1/Box
C7204	Probes with irrigation holes	10mm	33cm	1/Box
C7205	Pool tip probe with irrigation holes	5mm	33cm	1/Box

* For use with reusable probe

** For use with reusable probe and reusable tubing

Visit www.appliedmedical.com/epix
for more information

© 2018 Applied Medical Resources Corporation. All rights reserved.

Applied Medical, the Applied Medical logo design and marks designated with a ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, and/or the European Union. SC01715F

