



Latis[®]
Graft Cleaning System

latis

Graft Cleaning System

The Advanced Balloon Catheter

The *latis* graft cleaning catheter is the most advanced balloon catheter designed for removing clots in non-native vessels

The unique traction balloon engages tough adherent thrombi while maintaining tactile feedback for effective treatment


Effective Clot Removal

Latex-free symmetrical balloon maintains a cylindrical shape for improved clot management

Unique *latis* mesh allows for efficient cleaning and effective removal of clots from synthetic grafts

Precise Atraumatic Navigation

Flexible tip and traction balloon offer tactile feedback





latis Graft Cleaning System

Reorder No.	Description	Size	Length
A4G00	<i>latis</i> Graft Cleaning Catheter	4F	80cm
A4G02	<i>latis</i> Graft Cleaning Catheter	4F	40cm
A4GW6	<i>latis</i> over-the-wire	5F	60cm

Visit www.appliedmedical.com/latis
for more information