

STERILE SPECIMEN COLLECTION SWAB

6-Inch Honey Dipper Nasopharyngeal Swab
with 70mm Breakpoint



As we navigate through these unprecedented times, Applied Medical is committed to rapidly responding to the needs of the healthcare community. We are all in this together, fighting against COVID-19.

6-INCH HONEY DIPPER NASOPHARYNGEAL SWAB WITH 70mm BREAKPOINT

PRODUCT DESCRIPTION

Applied Medical's honey dipper sterile specimen collection swab is indicated for the collection and transportation of nasopharyngeal clinical samples to support COVID-19 testing. The swab consists of a medical-grade polystyrene handle and honey dipper tip.

FEATURES

- Sterile, individually packaged swabs
- Made in the USA with materials sourced domestically

DIMENSIONS

Overall Length: 5.9 inches

Tip

- Design: Honey dipper
- Material: Polystyrene
- Length: 0.715 inches
- Diameter: 0.014 inches

Handle

- Material: Polystyrene
- Length: 3.097 inches
- Diameter: 0.098 inches

Breakpoint: 70mm

ORDER INFORMATION

Model #	Description	Qty/Box
NPS19	Sterile specimen collection swab	200*

*Controlled by weight

Applied Medical is pleased to offer nasopharyngeal swabs direct to U.S. healthcare facilities and government agencies. For more information or to inquire about placing an order, please contact COVIDResponse@appliedmedical.com.

